



Advanced Chemistry for Advanced Practices



DARBY DENTAL SUPPLY, LLC



Opti-Cide³

Fungicidal • Bactericidal • Virucidal • Tuberculocidal

A 3-Minute Ready-To-Use Disinfectant, Decontaminant & Cleaner

- Equipment & Surface Friendly
- Fragrance-Free
- Dye-Free
- Blue Surgical-Grade Wipe
- Broad spectrum kill in three minutes

...That's the Opti-Cide³ Advantage!

Compatibility of Opti-Cide³ on Common Dental Surfaces

	Opti-Cide ³	Cavicide [®]	Discide [®] Ultra	Lysol [®] I.C.	Bleach
Smooth Vinyl	✓	✓	✓	✓	
Textured Vinyl	✓	✓	✓	✓	✓
Brushed Aluminum	✓		✓	✓	
Stainless Steel	✓				✓
Formica	✓				
Hosing	✓		✓	✓	
Plexiglas	✓	✓			
High Impact Plastic	✓				
Enameled Metal	✓				

✓ = No significant damage to surface after long-term repeated disinfection as measured by cleanliness – more complete removal of blood.

Data on file, Palenik, C.J., et. al., 2008, Effect of Long-term Disinfection on Dental Clinical Contact Surfaces.

Please email sciencegeeks@biotrol.com for complete abstract and results. Discide is a registered trademark of Palmero Healthcare. Cavicide is a registered trademark of Metrex Research Corporation[®]. Lysol is a registered trademark of Reckitt Benckiser Group.



BIOTROL™ Infection Control Down to a Science



To Reorder Call
800.645.2310 or
Visit darby.com



Opti·Cide³



Fungicidal • Bactericidal • Virucidal • Tuberculocidal

A 3-Minute Ready-To-Use Disinfectant, Decontaminant & Cleaner

DARBY DENTAL SUPPLY, LLC

Advanced Chemistry

Opti·Cide³ is a synergistic formula containing surfactants and dual quaternary ammonium actives.

Dual quats offer several advantages over previous generation single chain quats*:

- Improved biocidal activity
- Stronger detergency
- Relatively lower level toxicity than previous quat formulations

*Block, S.S. Disinfection, Sterilization, and Preservation. Lippincott, Williams and Wilkins: Philadelphia, 2001. 5th Ed; p 292.

Surgical Grade Wipes



- Heavily saturated
- Blue, 7 layer wipe
- Surgical-grade



Surgical-grade Opti·Cide³ wipes won't lint or leave residue, protecting your patients, personnel, and practice.

Product Packaging

PRODUCT	ORDER #
Solution:	
(EPA REG. NO. 70144-1-51003)	
24 oz. trigger spray bottle	8126204 (DOCS12-024)
1 gallon pour bottle	8126206 (DOCP04-128)
2.5 gallon with spigot	8126208 (DOCP02-320)

Surface Wipes:	
(EPA REG. NO. 70144-2-51003)	
100 each large (7" x 10") Wipes saturated with Opti·Cide ³ solution packed in a round dispenser canister	8126210 (DOCW06-100)



Description & Use

Opti·Cide³ is a ready-to-use cleaner that provides complete disinfection* and is specially formulated for use with all dental facilities and surfaces.

- Opti·Cide³ effectively and quickly loosens and suspends blood and bodily fluids, organic and non-organic soil including spores and other cellular debris for easy removal during the cleaning process.
- Opti·Cide³ is a pre-cleaning instrument spray. Place instruments into suitable container. Thoroughly spray instrument surfaces with undiluted Opti·Cide³ to thoroughly drench surfaces. Allow instruments to remain wet for 3 minutes. Thoroughly rinse, then follow with the appropriate terminal cleaning and sterilization/disinfection process.
- Opti·Cide³ may be used to treat dental impressions. Simply spray impression and allow to be wet for 3 minutes. Follow with a cool tap water rinse and drying.

*Opti·Cide³ effective against TB in 3 minutes. Please see label for full instructions.



BIOTROL™

Infection Control Down to a Science



To Reorder Call
800.645.2310 or
Visit darby.com