


# Hand Peace.

*With personality...*



DARBY DENTAL SUPPLY, LLC

PerfectPearl™ Brush, Classic & Midnight Series



*PerfectPearl Brush*  
DISPOSABLE PROPHY BRUSH

*PerfectPearl Classic*  
DISPOSABLE PROPHY ANGLE

*PerfectPearl Midnight*  
DISPOSABLE PROPHY ANGLE

**PerfectPearl™ disposable prophy angle – experience comfort, control and innovation in disposable prophy angles.**

## The PerfectPearl disposable prophy angle.

With its unique beveled gearing, these contra angles provide smooth, low vibration operation. The soft prophy cup offers nice flare for excellent subgingival cleaning and the spiral grooves help push paste towards the tooth, reducing splatter. Ask your dealer for **PerfectPearl** and start to experience **Perfection** today.

- **Low Vibration and Chatter.** Designed for all day comfort.
- **Exceptional interproximal access.** Easy stain removal.
- **Reduces splatter.** Spiral grooves on cup help push paste towards the tooth.
- **Super control.** Smooth beveled gears offer low vibration.
- **Latex free soft cup.** PerfectPearl does not contain natural rubber latex.
- **Soft grip pad.** Easy to hold and grip.
- **Economical.** Competitively priced.
- **Contra angle design.** Ergonomically designed for clinician comfort.
- **Small size in patient's mouth.** Good visibility speeds the procedure.



**KerrTotalCare**

*PerfectPearl*<sup>™</sup>  
 Disposable Prophy Angle



The contra angle design and straight shaft offer comfort and control, with low vibration.

Latex Free:  
 products do not contain natural or synthetic rubber latex.



Straight shaft designed to reduce chatter and vibration during procedure.

## PerfectPearl Order Information (US & CA)

### PerfectPearl™ Classic

Item #	Product Description	Packaging
9551595 (10-9510)	Dark Blue Pad, Standard Cup	200ct per box
9551593 (10-9511)	Dark Blue Pad, Standard Cup, Bulk Box	1200ct per box
8546272 (10-9512)	Pink Pad, Standard Cup	200ct per box
8546274 (10-9513)	Purple Pad, Standard Cup	200ct per box
8546276 (10-9514)	Dark Blue Pad, Soft Cup	200ct per box
8546278 (10-9515)	Pink Pad, Soft Cup	200ct per box
8546280 (10-9516)	Purple Pad, Soft Cup	200ct per box

### PerfectPearl™ Midnight

Item #	Product Description	Packaging
8546284 (10-9517)	Purple Pad, Standard Cup	200ct per box
8546286 (10-9518)	Pink Pad, Standard Cup	200ct per box
8546288 (10-9519)	Teal Pad, Standard Cup	200ct per box
8546290 (10-9521)	Purple Pad, Soft Cup	200ct per box
8546292 (10-9522)	Pink Pad, Soft Cup	200ct per box
8546294 (10-9523)	Teal Pad, Soft Cup	200ct per box

### PerfectPearl™ Brush

Item #	Product Description	Packaging
8546296 (10-9524)	Tapered Brush	200ct per box



To Reorder Call  
**800.645.2310**  
 or Visit [darby.com](http://darby.com)

**KerrTotalCare**