




DARBY DENTAL SUPPLY, LLC

ultrammat 2

HIGH SPEED AMALGAMATOR

HOCHFREQUENZ MISCHGERÄT
AMALGAMADOR DE ALTA VELOCIDADE
AMALGAMADOR DE ALTA VELOCIDAD
VIBREUR ELECTRONIQUE A GRANDE VITESSE
VIBRATORE PER AMALGAMA AD ALTA VELOCITÀ
HIGH SPEED CAPSULESCHUDDER
HØJHASTIGHEDS-KAPSELRYSTER
HØYHASTIGHETS-KAPSELRYSTER
NOPEA KAPSELISEKOITIN
ΥΨΗΛΗΣ ΤΑΧΥΤΗΤΑΣ ΔΟΝΗΤΗΣ ΑΜΑΛΓΑΜΑΤΟΣ
WYSOKOOBROTOWY MIESZALNIK
NAGY SEBESSÉGŰ KAPSZULARÁZÓGÉP
SUURE KIIRUSEGA AMALGAAMISEGISTI
HITRI MEŠALEC AMALGAMA
AUGSTAS FREKVENCES AMALGAMATORS
DIDELIŲ APSUKŲ AMALGAMŲ MAIŠYTUVAS
VYSOKORYCHLOSTNÍ AMALGAMÁTOR
HITRI MEŠALEC AMALGAMA
高速アマルガムミキサー
高速銀汞調拌机



SDI


ultramat 2

5 year warranty
amalgamator



- Safety feature
- Consistent accurate trituration
- Practical
- Multi-use mixer
- High energy output
- Low energy usage
- 5 year warranty
- Simple to use
- Quiet low vibration
- Dual voltage


5 years warranty



ultrammat 2 features:


Consistent and accurate trituration

The Ultramat 2 features a 0-16 second microprocessor-crystal controlled oscillator timer.




Multi-use mixer

The multi-use Ultramat 2 is suitable for all types of encapsulated dental materials.




Quiet low vibration

The Ultramat 2 is the quietest, low vibration mixer.




Safety feature

The micro-switch automatically halts trituration when the lid of the Ultramat 2 is lifted.




Simple to use

The touch control system of the Ultramat 2 is quick and simple to operate.




High energy output

The high energy Ultramat 2 has a mixing frequency of 4600 oscillations per minute +/-1%.




Practical

The one surface cover facilitates infection control cleaning and the non-porous plastic surface can be cleaned with any dental disinfectant.



Low energy usage

The Ultramat 2 has economical low energy usage.




5 year warranty

The Ultramat 2's robust and reliable construction assures your purchase is a long lasting investment.



Dual voltage

The Ultramat 2 operates on voltages between 120 and 240v, 50 and 60Hz, for compatibility with multiple power voltages and frequencies.



ultrammat 2



Dimensions:

width	depth	height
217mm	173mm	180mm
8 1/2"	6 3/4"	7 1/8"

Voltage:

120 / 240 volts
50 / 60 Hz

Weight:

2.765kg
6.095lbs

Ultramat 2

(USA - includes power cord for USA)

4473100 (5546045)

*AutoMix, Dentomat Compact, OptiMix, ProMix, RotoMix, Silamat S5, Wig-L-Bug, are not registered trademarks of SDI Limited.



Made in Australia by SDI Limited
 Bayswater, Victoria 3153
 Australia 1 800 337 003
 Austria 00800 022 55 734
 Brazil 0800 770 1735
 France 00800 022 55 734
 Germany 0800 100 5759
 Ireland 01 886 9577
 Italy 800 780625
 New Zealand 0800 734 034
 Spain 00800 022 55 734
 United Kingdom 00800 022 55 734
 USA & Canada 1 800 228 5166
www.sdi.com.au



DARBY DENTAL SUPPLY, LLC

To Reorder Call 800.645.2310
 or Visit darby.com



DENTAL EQUIPMENT



M554601 D

+DO36M55460112

ACN 008 075 581