

MICROFLEX®

▲ COBALT®



BLUE
NITRILE
NON-STERILE
AMBIDEXTROUS
FULLY TEXTURED
BEADED CUFF
FOR SINGLE USE ONLY
NON-LATEX

POWDER-FREE EXAMINATION GLOVES

This soft and stretchy nitrile glove is the best-selling glove from Microflex® brand of products. Durable and puncture resistant with a fully textured surface for enhanced grip. Ideal for jobs requiring exceptional tactile sensitivity.

- Unique, cobalt blue coloring
- Consistent product with exceptional fit and feel



PRODUCT CODE			PACKAGING	
3173220	XS	(N190)	100	Gloves/Box
3177321	S	(N191)		
3177322	M	(N192)		
3177323	L	(N193)		
3177324	XL	(N194)		

GLOVE PROFILE		ASTM D 6319	COBALT®
Length	<i>mm (inches)</i>	min. 220 (8.7)	240 (9.5)
Cuff Thickness	<i>mm (mils)</i>	n/a	0.09 (3.5)
Palm Thickness	<i>mm (mils)</i>	min. 0.05 (2.0)	0.10 (3.9)
Finger Thickness	<i>mm (mils)</i>	min. 0.05 (2.0)	0.12 (4.7)
Tensile Strength	<i>Before Aging</i>	min. 14	18
	<i>After Aging</i>	min. 14	14
Elasticity/Elongation	<i>Before Aging</i>	min. 500	500
	<i>After Aging</i>	min. 400	400
Pinhole Defect Rate*	<i>(AQL)</i>	max. 2.5	1.5

* Pinhole Defect Rate is FDA minimum requirement.



RECOMMENDED INDUSTRIES



DARBY DENTAL SUPPLY, LLC

To Reorder Call 800.645.2310
or Visit darby.com

Ansell

® and ™ are trademarks owned by Ansell Limited or one of its affiliates. © 2014 All Rights Reserved.
COLL-SPEC-N19. v1