




DARBY DENTAL SUPPLY, LLC

Sani-Tip Type Disposable Air/Water Syringe Tips



Darby's Disposable Air/Water Sani-Tip Type Tips are made of durable plastic. These tips feature a central water channel surrounded by six individual air channels, this design delivers instant switch between air and water resulting in consistent dry air.

Features and Benefits:

- **Clear Plastic Disposable Tip**
- **Eliminates Cross Contamination**
- **Provides Optimal Bonding Performance**
- **Decreases Chance Of Adhesive Failures**
- **Cost Effective**
- **Eliminates Moisture**
- **Fits In Current Sani-Tip Adaptor**

250/Bag

9528894



Darby Dental Supply, LLC

800-645-2310

www.darby.com