



DESIGN



TRUST



CONFIDENCE

Introducing **SECURE FIT**® Technology
THE DENTAL MASK OF CHOICE.



- ➔ Aluminum strips both above the nose and below the chin for a custom fit.
- ➔ Patent pending design significantly reduces gapping at top, bottom and sides.

Better Fit.

The only mask that contours to any size or shaped face for a custom fit.

Better Filtration.

Filters harmful bacteria, viruses and aerosols.

Better Protection.

3x more protective than standard face masks against airborne particulates.*

*Study on file



Top Mask



PROTECTING YOU AND YOUR PATIENTS.



Products manufactured by Crosstex

FACE MASKS WITH SECURE FIT[®] TECHNOLOGY

Better Fit. Better Filtration. Better Protection.

When compared to standard masks, Crosstex[®] face masks with Secure Fit[®] Technology are proven to provide up to **3 times** more filtration, fluid and aerosol protection!* They are available in 3 levels, so you always have the appropriate level of mask protection for the procedure you are performing.

*Study on File



3x's More Protection

- A patent-pending design reduces gapping at top, bottom and sides by incorporating aluminum strips both above the nose and under the chin.
- Custom fits to any size or shaped face to reduce your exposure to airborne particulates.



Aluminum nose and chin pieces (green arrows) reduce gapping (purple arrows).

3 Levels of Protection

3
99.9%
160^{mm}Hg

Level 3 - Heavy Amounts of Fluid

The Ultra[®] Line

Ideal for procedures where heavy amounts of fluid, spray or aerosols are produced.

2
99.9%
120^{mm}Hg

Level 2 - Moderate to Light Amounts of Fluid

The Procedural Line

Ideal for Procedures where moderate to light amounts of fluid, spray or aerosols are produced.

1
98.4%
80^{mm}Hg

Level 1 - Low Amounts of Fluid

The Isofluid[®] Line

Ideal for procedures where low amounts of fluid, spray or aerosols are produced.

Description	Color	Quantity	SKU Number	Ref. #
The Ultra[®] Line: Exceeds ASTM Level 3				
	Ultra [®] Earloop w/ Secure Fit [®] Technology	Blue	50/Box	3413166 (GCFCXUSF)
	Ultra [®] Sensitive Earloop w/ Secure Fit [®] Technology*	White	50/Box	3413175 (GCFCXSSF)
	Ultra [®] Sensitive FogFree [®] Earloop w/ Secure Fit [®] Technology*	White	40/Box	3413178 (GCFCXSFSF)
	Ultra [®] Sensitive FogFree [®] Earloop w/ Shield w/Secure Fit [®] Technology*	White	25/Box	3413179 (GCPWSSF)
* Ultra Sensitive masks are designed for sensitive skin: chemical free, hypo allergenic, and free of inks and dyes.				
The Procedural Line: Exceeds ASTM Level 2				
	Procedural Earloop w/ Secure Fit [®] Technology	Blue	50/Box	3413168 (GCPBLSF)
		Lavender	50/Box	3413167 (GCPLVSF)
		Pink	50/Box	3413169 (GCPPKSF)
The Isofluid[®] Line: Exceeds ASTM Level 1				
	Isofluid [®] Earloop w/ Secure Fit [®] Technology	Blue	50/Box	3413162 (GCIBLSF)
		Lavender	50/Box	3413163 (GCILVSF)
		Pink	50/Box	3413164 (GCIPKSF)
	Isofluid FogFree [®] Earloop w/ Secure Fit [®] Technology	Blue	40/Box	3413176 (GCICXBSF)
	Isofluid FogFree [®] Earloop w/Shield w/Secure Fit [®] Technology	Blue	25/Box	3413177 (GCIPWBSF)
NOTE: Filtration referenced is BFE (Bacterial Filtration Efficiency)				
Optically clear, distortion free, wrap-around Face Shield	Clear	24/Box	3411400	(GCSS)



To Reorder Call 800.645.2310
or Visit darby.com



PROTECTS[®]